

Akriti[®]

Medical Disposable Products



PRODUCTS RANGE

PAGES

SURGICAL GOWNS RANGE

1 - 2

PROTECTIVE APPARELS

3

FACE MASK & ACCESSORIES


4

SURGICAL DRAPES & PACKS

5 - 15

EQUIPMENT COVERS

16



Surgical Gowns

“Akriti” gowns has designed with safety, comfort and flexibility. We recognise that user may have different preferences for Surgical Gown. “Akriti”

offers Four Styles Each of Style offers different fabric that has its own unique, qualities.

Maxcare : SMS Fabric Gowns

Micromax : SMMS Fabric Gowns

Ultra care : ACT Fabric Gowns

Swiftcare: Dupont spunlace Gowns

Swiftsoft : BVB Fabric Gowns

Meets All International Standards

Akriti gowns meets and/or exceed the latest Meets All International Standards from AAMI (Association for the advancement of medical instrumentation) and EN13795 standard for surgical gowns.

Lowlinting

Surgical gowns are very low linting, minimising the risk of introducing particles into the O.R.

–Features of Akriti Gowns


- Neck – Adjustable with Velcro
- Medical Cuffs – Made from cotton.
- Latex free. No itching or irritation
- Breathable
- Impervious Fluid repellent
- Size – S / M / L / XL / XXL Conveniently folded
- Ultrasonic sealing for –maxcare, micromax,
- Ultracare and Medimax
- Sterile by ETO



Maxcare

SMS Fabric Gowns

Ahlstrom SMS fabrics are available with antistatic and alcohol-repellent finishes. They deliver high fluid and particulate barrier and are soft, strong, and low in lint. Compatible with ethylene oxide, steam, and gas plasma sterilization, Ahlstrom SMS fabrics are available for AAMI PB70 Level 3 and EN 13795 High Performance, less critical area and Standard Performance, all areas.




Standard Surgical Gown
(Sterile, Packed with 2 O.T. Towels)

Product Code	Size	Specification
MCS : 1048	Small	Standard
MCS : 1050	Medium	Standard
MCS : 1052	Large	Standard
MCS : 1054	Xtra Large	Standard
MCS : 1058	X-Xtra Large	Standard


Surgical Reinforced Gown
(Sterile, Packed with 2 O.T. Towels)

Product Code	Size	Specification
MCR : 1048	Small	Reinforced
MCR : 1050	Medium	Reinforced
MCR : 1052	Large	Reinforced
MCR : 1054	Xtra Large	Reinforced
MCR : 1058	X-Xtra Large	Reinforced





Wraparound Reinforced Gown
(Sterile, Packed with 2 O.T. Towels)

Product Code	Size	Specification
MCP : 1048	Small	PE Reinforced
MCP : 1050	Medium	PE Reinforced
MCP : 1052	Large	PE Reinforced
MCP : 1054	Xtra Large	PE Reinforced
MCP : 1058	X-Xtra Large	PE Reinforced



Protective Apparels

Scrub Suit




Product Code	Contents	Size
MCP 2010	Material - SMS/Spunlace	L
MCP 2011	Material - SMS/Spunlace	XL

Coverall Suit

Product Code	Contents	Size
MCP 2020	Material - SMS/Spunbond	XL
MCP 2021	Material - SMS/Spunbond	XXL

Lab Coat

Product Code	Contents	Size
MCP 2030	Material - SMS/Spunbond	Medium
MCP 2031	Material - SMS/Spunbond	L
MCP 2032	Material - SMS/Spunbond	XL




Patient Gown

Product Code	Contents	Size
MCP 2040	Material - SMS/Spunbond	Small
MCP 2041	Material - SMS/Spunbond	Medium
MCP 2042	Material - SMS/Spunbond	L
MCP 2043	Material - SMS/Spunbond	XL

Protective Apron

Product Code	Contents	Size
MCP 2050	Material - Plastic	Small
MCP 2051	Material - Non-Woven	Medium
MCP 2052	Material - PVC	L

Protective Apron



Product Code	Contents	Size
MCP 2060	Material - Spunbond	Small
MCP 2061	Material - Spunbond	Medium
MCP 2062	Material - Spunbond	L

Face Mask & Accessories



Mediblu Face Masks Features

- Made of Spunbond pp fabric
- Mono / Bi / Three layered
- Made with International Ultrasonic Sealing Technology
- Option of Tie & Elastic



Type Of Face Masks

- 1 Ply Face Mask - Tie/earloop
- 2 Ply Face Mask - Tie/earloop
- 3 Ply Face Mask - Tie/earloop
- Packing Per Box 100 Pieces



Mediblu Surgeon Caps

- Surgeon Caps - Tie/Elastic
- Bouffant Cap
- Nurses Cap
- Shoe Cover - Plastic/Non Woven



Other Accessories

- Eye Shield Face Mask
- Disposable Goggles
- Hood With Eye Shield
- Shoe Cover - Plastic/Non Woven
- Disposable Hood Cap
- Camera Cover
- Disposable Panties - Spunlace/Spunbond
- Disposable Facial Gown



Surgical drapes are designed to establish and maintain a sterile field. It's important to select a surgical drape that provides reliable performance. To achieve the best surgical outcomes look for surgical drapes with qualities like impervious to strike through, fluid control feature, reliable adhesion, superior strength and comfort.

Design

Mediblu drapes are designed, folded and packed as recommended by The Association of Perioperative Registered Nurses (AORN). Clearly labeled instructions on all drapes to facilitate correct, aseptic and fast application. Wider anaesthesia screen and arm board cover to drape a person from fingertip to fingertip.

Reinforcement

Large impervious reinforcement zone surrounded with fenestrations & splits of each drape. Our reinforcement fabric is highly absorbent to absorb fluids and blood and is plastic backed.

Fluid Containment

In built pouches offered in selected drapes to collect & contain blood and body fluids, minimizing the risk of exposure. Pouches also allow for efficient, easy clean up & fluid management.

Adhesive & Incise

We use premium quality, skin friendly adhesive tape to secure drapes in place without compromising barrier integrity. Incise is included in many procedural drapes such as C- Section drapes, craniotomy drape, open heart drape etc.

Tube & Line Holders

Mediblu drapes feature Velcro style tube and line holders to place surgical lines and tubes and help keep them clear of the fenestration.

Surgical Drapes and Packs



Highlights of Mediblue Drapes

- ▲ Low-linting, non woven, non-absorbent fabric
- ▲ Absorbent only around the fenestration, surgical and critical zone areas, if applicable
- ▲ Soaks up spills and splatters
- ▲ No fluid strike-through
- ▲ Made from absorbent and repellent materials like Kraton, SMS, Polypropylene
- ▲ Continuous broad-spectrum antimicrobial activity in the drape adhesive where it can't be washed away
- ▲ Clinically proven to help reduce the risk of wound contamination and bacterial migration
- ▲ No pinning or clips required

Choice Of Fabrics

In Mediblue you have a choice of surgical drapes and packs in different fabrics available as SMS, Spun lace, laminate without any compromise in quality.

Our range of products include disposable surgical drapes, Disposable surgical Packs, Disposable drapes sheet, sterile surgical drape and many more.

SMMS non-woven fabric

Outerlayer


Spunbond repellent fabric

2 Middle layers

Meltblown, fluid control, bacterial barrier

Innerlayer

Spunbond strength




- Lightweight for maximum comfort
- High barrier performance
- Good strength characteristics
- Good breathability
- Low linting




Surgical Packs (Orthopedic)

Product Code : MB - TKR - M - 01 **Total Knee Replacement Pack**




Quantity	Contents	Size
1 Pc.	Main "O" Drape	345 x 180 x 250
1 Pcs.	Top Sheet	200 x 150 cms.
1 Pc.	Bottom Drape	275 x 180 cms.
1 Pc.	Trolley Sheet	210 x 163 cms.
4 Pcs.	Wr. around Gown Reinforced	XL

Product Code : MB - TKR - M - 02 **Total Knee Replacement Pack**




Quantity	Contents	Size
1 Pc.	Main "O" Drape	345 x 180 x 250
1 Pcs.	Top Sheet	200 x 150 cms.
1 Pc.	Bed Sheet	160 x 225 cms.
1 Pc.	Disposable Bag	30 x 40 cms.

Product Code : MB - THR - M - 01 **Total Hip Replacement Pack**




Quantity	Contents	Size
1 Pc.	Main "U" Drape	300 x 200 cms.
1 Pcs.	Top Sheet	200 x 150 cms.
1 Pc.	Bottom Drape	275 x 180 cms.
1 Pc.	Trolley Sheet	210 x 163 cms.
1 Pc.	Disposable Bag	210 x 163 cms.
4 Pcs.	Wr. around Gown Reinforced	XL

Product Code : MB - THR - M - 02 **Total Hip Replacement Pack**



Quantity	Contents	Size
1 Pc.	Main "U" Drape	340 x 196 x 260
1 Pcs.	Top Sheet	200 x 150 cms.
1 Pc.	Bed Sheet	160 x 225 cms.
1 Pc.	Disposable Bag	15 x 30 cms.







Surgical Packs (Cardiovascular)

Product Code : MB - CABGP - M01

CABG Pack 01


Quantity	Contents	Size
2 Pcs.	Main Drape	250 x 120 cms
1 Pcs.	Bottom Drape	200 x 175 cms.
4 Pcs.	Absorbent Towel	38 x 38 cms.
4 Pcs.	Square Sheet	100 x 100 cms.
1 Pc.	Anaesthesia Isolation Sheet	150 x 300 cms.
2 Pcs.	Waist Band	100 x 50 cms
1 Pair	Shoe Cover	38 x 18 cms.
2 Pcs.	Tubing Organizer	20 x 30 cms.
1 Pc.	Disposable Bag	30 x 45 cms.
1 Pc.	Mayo Trolley Cover	75 x 100 cms.
2 Pcs.	Trolley Cover	250 x 200 cms.
2 Pcs.	Basin Stand Cover	225 x 175 cms.
2 sets	Surgical Gown XL & L each	XL & L



Product Code : MB - CABGP - M02

CABG Pack 02


Quantity	Contents	Size
2 Pcs.	Large Sheet	250 x 160 cms
2 Pcs.	Draw Sheet With Adhesive Tape	200 x 100 cms.
1 Pc.	Draw Sheet	160 x 120 cms.
1 Pc.	Anaesthesia Isolation Sheet	235 x 160 cms.
1 Pc.	Square Sheet	80 x 50 cms
2 Pcs.	Waist Band	60 x 45 cms.
4 Pcs.	Small Towel	30 x 30 cms.
1 Pc.	Disposable Bag	40 x 30 cms.
1 Pc.	Main Drape With 3M Ioban	250 x 180 cms
1 Pc.	Mayo Stand Cover	150 x 60 cms.
2 Pcs.	Basin Stand Cover	150 x 100 cms.
1 Pc.	Trolley Cover	200 x 160 cms.
2 Pcs.	Big Towel	75 x 75 cms
2 Pcs.	Small Tubing Organizer	10 x 15 cms.
1 Pc.	Big Tubing Organizer	20 x 30 cms.
2 Pcs.	Wraparound Gown With Towel	XL
2 Pcs.	Wraparound Gown With Towel	L




Surgical Packs (Cardiovascular)

Product Code : MB - CABG - M03


CABG Pack III



PLD Sheet



Abdominal Drape



Mayo Cover



Surgeon Gown

Quantity	Contents	Size
1 Pc.	PLD/Plastic Sheet	125 x 150 cms.
1 Pcs.	PLD/Plastic Sheet	125 x 125 cms.
1 Pc.	Groin Skin Towel	65 x 65 cms.
1 Pc.	Abdominal Drape	86 x 36 cms.
3 Pcs.	Leg Skin Towel	87 x 87 cms.
1 Pc.	Leg Cover - Plain	140 x 140 cms.
1 Pc.	Leg Cover - with Release	140 x 140 cms.
2 Pcs.	Side Drapes	200 x 107 cms.
1 Pc.	Cardio Thorasic Drape	300 x 310 cms.
4 Pcs.	Cauterg Bag	28 x 33 cms.
1 Pc.	Pump Sheet	180 x 120 cms.
1 Pc.	Release Tapes	18 x 5 cms.
Mayo Cover Kit		
1 Pc.	Mayo Cover	150 x 60 cms.
1 Pc.	Instrument Holder	160 x 80 cms.
Surgeon Gown Kits		
4 Pcs.	Surgeon Gowns	Medium
2 Pcs.	Surgeon Gowns	XL
6 Pcs.	Wipes	35 x 50 cms.
1 Pc.	Wrap	100 x 100 cms.
Table Cover Kit		
2 Pcs.	Main Table Cover Large	310 x 180 cms.
2 Pcs.	Main Table Cover Medium	225 x 135 cms.
1 Pc.	Main Table Cover Small	90 x 90 cms.
1 Pc.	Instrument Holder	160 x 80 cms.
1 Pc.	Cardial Cover	28 x 60 cms.
1 Pc.	Suture Cover	65 x 65 cms.
1 Pc.	Wrap	80 x 80 cms.
Basin Stand Cover Kit		
2 Pcs.	Basin Stand	67 x 67 cms.
2 Pcs.	Basin Sheet	150 x 100 cms.
1 Pc.	Wrap	60 x 60 cms.
Redial Arm Pack		
1 Pc.	PLD	125 x 125 cms.
1 Pc.	Hand Skin Towel	65 x 65 cms.
1 Pc.	Hand Drape	81 x 81 cms.




Surgical Packs (Gynaecology)


Product Code : MB - GYP - M01

Gynaecology Pack


Quantity	Contents	Size
1 Pc.	Instrument Table Cover	140 x 160 cms
1 Pc.	Under Buttock Drape	35 x 100 cms.
1 Pair	Leggings	90 x 40 cms.
1 Pc.	T Shaped Drape	135 x 200 cms.
2 Pcs.	Adhesive Tape	8 x 45 cms.



Under Buttock Drape
with Fluid Collection Pouch



Abdominal Drape



Leggings

Product Code : MB - DLV - M01


Delivery Kit - I

Quantity	Contents	Size
1 Pc.	Setup Cover	110 x 195 cms
1 Pc.	Under Buttock Drape With Fluid Collection Pouch	75 x 125 cms.
1 Pc.	Abdominal Drape Absorbent Laminate Reinforcement	235 x 160 cms.
1 Pc.	Baby Blanket	60 x 90 cms.
2 Pcs.	Leggingsl 15 cms. Cuff	75 x 120 cms.
2 Pcs.	Absorbent Towel	40 x 40 cms.
10 Pcs.	Guaze Swab X-Ray Detectable	10 x 15 cms
1 Pc.	Cord Clamp	Standard
1 Pc.	Bulb Syringe	Standard
1 Pc.	Vaginal Sponge	Standard
1 Pc.	Maternity Pad	Standard
1 Pc.	Placenta Basin	Standard
1 Pc.	Black Table Cover Reinforced	140 x 200 cms.






Product Code : MB - DLV - M02

Delivery Kit - II

Quantity	Contents	Size
1 Pc.	Setup Cover	110 x 195 cms
1 Pc.	Under Buttock Drape With Fluid Collection Pouch	75 x 125 cms.
1 Pc.	Abdominal Drape Absorbent Laminate Reinforcement	100 x 150 cms.
2 Pcs.	Leggingsl 15 cms. Cuff	75 x 120 cms.
1 Pc.	Black Table Cover Reinforced	140 x 200 cms.



Surgical Packs (General)



Product Code : MB - UP - M01

Universal Pack I

Quantity	Contents	Size
1 Pc.	Instrument Table Cover	140 x 190 cms
1 Pcs.	Mayo Stand Cover	140 x 70 cms.
4 Pcs.	Side Drape	115 x 100 cms.
1 Pc.	Top Drape	160 x 275 cms.
1 Pc.	Bottom Drape	225 x 160 cms.
6 Pcs.	Absorbent Towel	35 x 35 cms

Product Code : MB - UP - M02


Universal Pack II

1 Pc.	Mayo Stand Cover Reinforced	60 x 140 cms
1 Pc.	Top Drape	130 x 225 cms.
1 Pc.	Draw Sheet	160 x 120 cms.
	<ul style="list-style-type: none"> • Absorbent Laminate Reinforcement • Adhesive Area • Tube Holders 	
2 Pcs.	Side Drape	75 x 90 cms.
	<ul style="list-style-type: none"> • Absorbent Laminate Reinforcement • Adhesive Area 	
1 Pc.	Bottom Drape	150 x 180 cms
	<ul style="list-style-type: none"> • Absorbent Laminate Reinforcement • Adhesive Area • Tube Holders 	
1 Pc.	Suture Bag	20 x 30 cms.
1 Pc.	Back Table Cover Reinforced	140 x 200 cms.

Product Code : MB - UP - M03


Universal Pack III

Quantity	Contents	Size
1 Pc.	Mayo Stand Cover Reinforced	60 x 140 cms
4 Pcs.	Utility Drape With Adhesive	140 x 70 cms.
1 Pc.	Suture Bag	40 x 65 cms.
1 Pc.	Back Table Cover Reinforced	20 x 30 cms.




Surgical Drapes


Total Knee Replacement Drape




Total Knee Replacement Drape



Total Hip Replacement Drape



Leg O Drape



Leg Drape

Product Code

Contents

Size

MB - TKRD - M1	Total Knee Replacement Drape	325 x 185 x 260 cms
	<ul style="list-style-type: none"> • 10 cms. Conformable Fenestration • Absorbent Laminate Reinforcement • Arm Board Cover • Tube Holders 	

Total Hip Replacement Drape

Product Code

Contents

Size

MB - THRD - M1	Total Hip Replacement Drape	325 x 185 x 260 cms
	<ul style="list-style-type: none"> • 100 x 10 cms. Adhesive Split • Absorbent Laminate Reinforcement • Arm Board Cover • Tube Holders 	

Leg O Drape

Product Code

Contents

Size

MB - LOD - M1	Leg O Drape	250 x 180 cms
	<ul style="list-style-type: none"> • 10 cms. Conformable Fenestration • Absorbent Laminate Reinforcement • Tube Holders 	


Leg U Drape

Product Code

Contents


Size

MB - LUD - M1	Leg U Drape	250 x 180 cms
	<ul style="list-style-type: none"> • 100 x 10 cms. Adhesive Split • Absorbent Laminate Reinforcement • Tube Holders 	



Surgical Drapes


Angiography Drape



Angiography Drape

Product Code	Contents	Size
MB - ANGD - M1	Angiography Drape <ul style="list-style-type: none"> • 10 cms. 2 Conformable Fenestration • Absorbent Laminate Reinforcement • Tube Holders 	310 x 180 cms


Bilateral Drape



Bilateral Drape

Product Code	Contents	Size
MB - BLD - M1	Bilateral Drape <ul style="list-style-type: none"> • 10 cms. 2 Conformable Fenestration • Absorbent Laminate Reinforcement • Tube Holders 	300 x 210 cms


Arthroscopy Drape



Arthroscopy Drape

Product Code	Contents	Size
MB - ARD - M1	Arthroscopy Drape <ul style="list-style-type: none"> • 10 cms. Conformable Fenestration • Absorbent Laminate Reinforcement • Tube Holders 	225 x 190 cms

Spinal Art Drape




Spinal Art Drape

Product Code	Contents	Size
MB - SAD - M1	Spinal Art Drape <ul style="list-style-type: none"> SAD I - 25 x 60 cms. Split Adhesive SAD II - 25 x 25 cms. Split Adhesive SAD III - 15 x 25 cms. Split Adhesive SAD IV - 15 x 25 cms. Split Adhesive • Absorbent Laminate Reinforcement • Tube Holders 	300 x 210 cms



Surgical Drapes


Laproscopy T - Shape Drape



Laproscopy T- shaped Drape

Product Code	Contents	Size
MB - LSD - M1	Laproscopy T - Shape Drape <ul style="list-style-type: none"> • 20 x 15 cms. Split Adhesive • Absorbent Laminate Reinforcement • Tube Holders 	275 x 85 x 175 cms


CABG Drape



CABG Drape

Product Code	Contents	Size
MB -CABGD - M1	CABG Drape <ul style="list-style-type: none"> • 30 x 40 cms. with 3M Ioban • Absorbent Laminate Reinforcement • Tube Holders 	250 x 180 cms


Craniotomy Drape



Craniotomy Drape

Product Code	Contents	Size
MB -CRD - M1	Craniotomy Drape <ul style="list-style-type: none"> • 30 x 20 cms. Split Adhesive • Absorbent Laminate Reinforcement • Tube Holders 	250 x 180 cms


Gynecology Drape



Gynecology Drape

Product Code	Contents	Size
MB -GYND - M1	Gynecology Drape <ul style="list-style-type: none"> • 10 x 23 cms. Split Adhesive • Absorbent Laminate Reinforcement • Tube Holders 	210 x 80 x 135 cms

Ophthalmic Drape



Ophthalmic Drape

Product Code	Contents	Size
MB -OPD - M1	Ophthalmic Drape <ul style="list-style-type: none"> • 70 x 90 cms. • 50 x 50 cms. • 100 x 100 cms. • 10 x 10 cms. Split With In Size • 15 x 20 cms. Collection Pouch • Stiffner 	



Head Office
Akriti Ophthalmic LLC,
330 N BRAND BLVD STE 700
GLENDALE, CA 91203
USA
+1 (864) 624 6898
sales@akritiophthalmic.com
www.akritiophthalmic.com

Asia
Akriti Ophthalmic Pvt Ltd,
M34, M35, M36 Medical Devices Park,
Sultanpur Pathancheru,
Hyderabad - 502319,
Telangana State, India.
+91 995 954 1995
sales@akriti.co.in
www.akriti.co.in

Middle East & Africa
Akriti Medical Equipments Trading Co L.L.C
The Prism Tower,
Flat No 3201, Business bay, Dubai,
United Arab Emirates
+971 547 917 448
sales@akritigroup.com
www.akritigroup.com